

Product datasheet for **RC208216L4V**

SLITRK1 (NM_052910) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLITRK1 (NM_052910) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLITRK1
Synonyms:	LRRC12; TTM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_052910
ORF Size:	2088 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208216).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_052910.1
RefSeq Size:	5201 bp
RefSeq ORF:	2091 bp
Locus ID:	114798
UniProt ID:	Q96PX8
Cytogenetics:	13q31.1
Protein Families:	Transmembrane
MW:	77.6 kDa


[View online »](#)

Gene Summary:

This gene encodes a member of the SLITRK protein family. Members of this family are integral membrane proteins that are characterized by two N-terminal leucine-rich repeat (LRR) domains and a C-terminal region that shares homology with trk neurotrophin receptors. However, the protein encoded by this gene lacks the region of homology to neurotrophin receptors. This protein is thought to be involved in neurite outgrowth. Mutations in this gene may be associated with Tourette syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]