

Product datasheet for **RC208220A1V**

Human RNF11 (NM_014372) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RNF11 (NM_014372) AAV Particle
Tag: Myc-DDK
Symbol: RNF11
Synonyms: CGI-123; SID1669
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208220 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGGAAGTGCCTCAAATCCCCACCTCGGATGACATCTCCCTGCTTCACGAGTCTCAGTCCGACCGGG
CTAGCTTTGGCGAGGGGACGGAGCCGGATCAGGAGCCGCCGCCATATCAGGAACAAGTTCAGTTCC
AGTCTACCACCAACACCTAGCCAGACTCGGCTAGCAACTCAGCTGACTGAAGAGGAACAAATTAGGATA
GCTCAAAGAATAGGTCTTATAACAACATCTGCCTAAAGGAGTTTATGACCCTGGAAGAGATGGATCAGAAA
AAAAGATCCGGGAGTGTGTGATCTGTATGATGGACTTTGTTTATGGGGACCAATTCGATTCTGCCGTG
CATGCATCTATCACCTGGACTGTATAGATGACTGGTTGATGAGATCCTTCACGTGCCCTCCTGCATG
GAGCCAGTTGATGCAGCACTGCTTTCATCCTATGAGACTAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208220 protein sequence
Red=Cloning site Green=Tags(s)

MGNCLKSPTSDDISLLHESQSDRASFGEGTEPDQEPPIPPYQEQVPVYHPTPSQTRLATQLTEEEQIRI
AQRIGLIQHLPGVYDPGRDGSEKKIRECVICMMDFVYGDPPIRFLPCMHIYHLCIDDWLMRSFTCPSCM
EPVDAALLSSYETN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_014372
ORF Size:	462 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014372.3</u>
RefSeq Size:	3082 bp
RefSeq ORF:	465 bp
Locus ID:	26994
UniProt ID:	<u>Q9Y3C5</u>
Cytogenetics:	1p32.3
MW:	17.4 kDa