

## Product datasheet for **RC208228L4V**

### HEATR7B2 (MROH2B) (NM\_173489) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HEATR7B2 (MROH2B) (NM_173489) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HEATR7B2
Synonyms:	HEATR7B2; SPIF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_173489
ORF Size:	4755 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208228).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173489.4</a>
RefSeq Size:	5307 bp
RefSeq ORF:	4758 bp
Locus ID:	133558
UniProt ID:	<a href="#">Q7Z745</a>
Cytogenetics:	5p13.1
MW:	180.8 kDa
Gene Summary:	May play a role in the process of sperm capacitation.[UniProtKB/Swiss-Prot Function]


[View online »](#)