

Product datasheet for **RC208237L2V**

SLC44A1 (NM_080546) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC44A1 (NM_080546) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC44A1
Synonyms:	CD92; CDW92; CHTL1; CONATOC; CTL1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_080546
ORF Size:	1971 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208237).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080546.3
RefSeq Size:	10521 bp
RefSeq ORF:	1974 bp
Locus ID:	23446
UniProt ID:	Q8WWI5
Cytogenetics:	9q31.1-q31.2
Domains:	DUF580
Protein Families:	Transmembrane



[View online »](#)

MW: 73.3 kDa

Gene Summary: Choline transporter. May be involved in membrane synthesis and myelin production.
[UniProtKB/Swiss-Prot Function]