

Product datasheet for RC208251L1

OLAH (NM_001039702) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OLAH (NM_001039702) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	OLAH
Synonyms:	AURA1; SAST; TE2; THEDC1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208251).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001039702
ORF Size:	795 bp




[View online >](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001039702.1
RefSeq Size:	1675 bp
RefSeq ORF:	798 bp
Locus ID:	55301
UniProt ID:	Q9NV23
Cytogenetics:	10p13
Protein Pathways:	Fatty acid biosynthesis, Metabolic pathways
MW:	29.9 kDa
Gene Summary:	Contributes to the release of free fatty acids from fatty acid synthase (FASN). Has broad substrate specificity, giving rise to a range of free fatty acids with chain lengths between 10 and 16 carbon atoms (C10 - C16).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC208251L1



Double digestion of RC208251L1 using SgfI and MluI