

Product datasheet for **RC208254L4V**

RNF144 (RNF144A) (NM_014746) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RNF144 (RNF144A) (NM_014746) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNF144
Synonyms:	hUIP4; RNF144; UBCE7IP4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014746
ORF Size:	876 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208254).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014746.2
RefSeq Size:	5755 bp
RefSeq ORF:	879 bp
Locus ID:	9781
UniProt ID:	P50876
Cytogenetics:	2p25.1
Domains:	IBR
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 32.7 kDa

Gene Summary: This gene encodes a member of a family of RING finger domain-containing E3 ubiquitin ligases that also includes parkin and parc. The expression of this gene is induced by DNA damage. The encoded protein interacts with the cytoplasmic DNA-dependent protein kinase, catalytic subunit (DNA-PKcs) and promotes its degradation through ubiquitination. The orthologous mouse protein has been shown to interact with a ubiquitin-conjugating enzyme involved in embryonic development. [provided by RefSeq, Mar 2017]