

## Product datasheet for **RC208259A1V**

### Human **BLCAP (NM\_006698) AAV Particle**

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human BLCAP (NM\_006698) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** BLCAP  
**Synonyms:** BC10  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC208259 representing NM\_006698  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGTATTGCCTCCAGTGGCTGCTGCCCGTCTCCTCATCCCCAAGCCCCTCAACCCCGCCCTGTGGTTCA  
 GCCACTCCATGTTTCATGGGCTTCTACCTGCTCAGCTTCCTCCTGGAACGGAAGCCTTGACAATTTGTGC  
 CTTGGTTTTCTGGCAGCCCTGTTCTTATCTGCTATAGCTGCTGGGGAAACTGTTTCCTGTACCACTGC  
 TCCGATCCCCGCTTCCAGAATCGGCGCATGATCCCGCGTTGTGGGCACC

**ACGCGT**ACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC208259 representing NM\_006698  
 Red=Cloning site Green=Tags(s)

MYCLQWLLPVLLIPKPLNPALWFSHSMFMGFYLLSFLLEKPKPTICALVFLAALFLICYSCWGNCFLYHC  
 SDSPLPESAHDPGVVGT

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_006698  
**ORF Size:** 261 bp  
**Buffer:** PBS with 0.001% Pluronic F68



[View online »](#)

<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_006698.2</a></u>
<b>RefSeq Size:</b>	2057 bp
<b>RefSeq ORF:</b>	264 bp
<b>Locus ID:</b>	10904
<b>UniProt ID:</b>	<u><a href="#">P62952</a></u>
<b>Cytogenetics:</b>	20q11.23
<b>MW:</b>	9.7 kDa