

Product datasheet for **RC208262L4V**

RAD21 (NM_006265) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RAD21 (NM_006265) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RAD21
Synonyms:	CDLS4; hHR21; HR21; HRAD21; MCD1; MGS; NXP1; SCC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006265
ORF Size:	1893 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208262).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006265.1
RefSeq Size:	3773 bp
RefSeq ORF:	1896 bp
Locus ID:	5885
UniProt ID:	O60216
Cytogenetics:	8q24.11
Domains:	Rad21_Rec8, Rad21_Rec8_N
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Cell cycle

MW: 71.7 kDa

Gene Summary: The protein encoded by this gene is highly similar to the gene product of *Schizosaccharomyces pombe rad21*, a gene involved in the repair of DNA double-strand breaks, as well as in chromatid cohesion during mitosis. This protein is a nuclear phospho-protein, which becomes hyperphosphorylated in cell cycle M phase. The highly regulated association of this protein with mitotic chromatin specifically at the centromere region suggests its role in sister chromatid cohesion in mitotic cells. [provided by RefSeq, Jul 2008]