

## Product datasheet for **RC208265L2V**

### **DRAK1 (STK17A) (NM\_004760) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	DRAK1 (STK17A) (NM_004760) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DRAK1
Synonyms:	DRAK1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004760
ORF Size:	1242 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208265).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004760.1</a> , <a href="#">NP_004751.1</a>
RefSeq Size:	2641 bp
RefSeq ORF:	1245 bp
Locus ID:	9263
UniProt ID:	<a href="#">Q9UEE5</a>
Cytogenetics:	7p13
Domains:	pkinese, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

**MW:** 46.4 kDa

**Gene Summary:** This gene is a member of the DAP kinase-related apoptosis-inducing protein kinase family and encodes an autophosphorylated nuclear protein with a protein kinase domain. The protein has apoptosis-inducing activity. [provided by RefSeq, Jul 2008]