

Product datasheet for **RC208271**

RAP1GAP (NM_002885) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RAP1GAP (NM_002885) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	RAP1GAP
Synonyms:	RAP1GA1; RAP1GAP1; RAP1GAPII; RAPGAP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC208271 ORF sequence

Red=Cloning site Blue=ORF Green=Tags(s)

 TTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGATTGAGAAGATGCAGGGAAGCAGGATGGATGAACAACGCTGCTCCTCCCGCCGCCCTCAAAACAG
 AGGAGGACTACATTCCATACCCGAGCGTGACGAGGTCTTGGGGCGAGAAGGACCCCTCCCCCTCATCCT
 GCTGCCCCAGTTTGGGGCTACTGGATTGAGGGCACCAACCACGAAATCACCAGCATCCCGAGACAGAG
 CCACTGCAGTCGCCCACAACCAAGGTGAAGCTCGAGTGCAACCCACAGCCCGCATCTACCGGAAGCACT
 TTCTCGGCAAGGAGCATTTCATTACTACTCACTGGACGCTGCCCTCGGCCACCTTGCTTCTCACTCAA
 GTACGATGTCATCGGGACCAAGAGCACCTGCGGCTGCTGCTCAGGACCAAGTGCCGGACATACCATGAT
 GTCATCCCATCTCCTGCCTACCGAGTTCCTAATGTTGTCCAGATGGCAAAGTTGGTGTGTGAAGACG
 TCAATGTGGATCGTTCTATCCTGTGCTCTACCCCAAGGCTTCCCGGCTCATCGTCACCTTTGACGAGCA
 TGTCATCAGCAATAACTTCAAGTTTGGCGTCATTTATCAGAAGCTTGGGCAGACCTCCGAGGAAGAACTC
 TTCAGCACCAATGAGGAAAGTCCCGCTTTCGTGGAGTTTCCTGAATTTCTTGGCCAGAAGGTCAAACCTGC
 AGGACTTTAAGGGGTTCGAGGAGGCCTGGACGTGACCCACGGGCAGACGGGGACCGAATCTGTGTACTG
 CAACCTCCGCAACAAGGAGATCATGTTTCACGTGTCCACCAAGCTGCCATACACGGAAGGGGACGCCAG
 CAGTTGCAGCGGAAGCGGCACATCGGGAACGACATCGTGGCTGTGGTCTTCCAGGATGAGAACACTCCTC
 TCGTGCCCGACATGATCGCGTCCAATTCCTGCATGCCTACGTCGTGGTGAGGCTGAGGGCGGGGGCCC
 TGATGGCCCCCTCTACAAGGTCTCTGCTCACTGCAAGAGATGATGTGCCCTCTTTGGACCCCCCTCCCG
 GACCCCGCTGTGTTCAAGGAAGGGCCTGAGTTCAGGAATTTTGTGACAAAGCTGATCAATGCTGAAT
 ATGCCTGCTACAAGGCAGAGAAGTTTGCCAAACTGGAGGAGCGGACGCGGGCCGCCCTCCTGGAGACGCT
 CTATGAGGAACTACACATCCACAGCCAGTCCATGATGGGCTTGGGCGGCGACGAGGACAAGATGGAGAAT
 GGCAGTGGGGGCGGCGGCTTCTTTGAGTCTTCAAGCGGGTCATCCGGAGCCGCAGCCAGTCCATGGATG
 CCATGGGGCTGAGCAACAAGAAGCCCAACACCGTGTCCACCAGCCACAGCGGGAGCTTCGCGCCCAACA
 CCCCAGCTGGCCAAGGCGGCTGGAATATCACTGATTGTCCCTGGGAAGAGCCCCACGAGGAAGAAGTCG
 GGCCCGTTTCGGCTCCCGCCGACGAGCGCCATTGGCATCGAGAACATACAGGAGGTGCAGGAGAAGAGGG
 AGAGCCCTCCGGCTGGTCAGAAGACCCAGACAGCGGGCACGTCTCACAGGAGCCCAAGTCGGAGAAGTC
 ATCCACTCAGAGCTCCCAGAGATGCCACGACCAAGAAGAGCGGAGACCGCAGCGCAGAGAGCAGAG
 GCGCTCAAGGACTTCTCCGCTCCTCGTCCAGTGCCAGCAGCTTCGCCAGCGTGGTGGAGGAGACGGAGG
 GTGTGGACGAGAGGACACAGGCCTGGAGAGCGTGTATCCTCAGGAACACCCCAAGCGGGACTCCTT
 CATCTATAGCACGTGGCTGGAGGACAGTGTGAGCACCAGTGTGGGGGACGCTCCCGAGGCCCTCTCGA
 TCACCCACCCAGACGCCGGAAGTTGGGGGACCTGCGTGTCCCGAGATCAAGATCCAGCTGGAAGCAT
 CTGAGCAGCACATGCCCCAGCTGGGCTGT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208271 protein sequence
 Red=Cloning site Green=Tags(s)


MIEKMQGSRMDEQRCSPFPLKTEEDYIPYPSVHEVLGREGPFPLILLPQFGGYWIEGTNHEITSIPETE
 PLQSPTTKVKLECNPTARIYRKHFLGKEHFNYYSLDAALGHLVFSKYDVIDGQEHRLLLRTKCRTHD
 VIPISCLTEFPNVVQMAKLVCEDEVNDRFYPLYPKASRLIVTFDEHVISNNFKFVVIYQKLQGTSEEL
 FSTNEESPAFVEFLEFLGQKVKLQDFKGRGLDVTHGQTGESVYCNRNKEIMFHVSTKLPTYEGDAQ
 QLQRKRHIGNDIVAVVFQDENTPLVPDMIASNFLHAYVVVQAEGGPDGPLYKVSVTARDDVPFFGPPLP
 DPAVFRKGPEFQEFLTLKLIAYACYKAEKFAKLEERTRAALLETLYEELHIHSQSMGLGGDEDKMEN
 GSGGGGFESFKRVIRSRSQSMDAMGLSNKKPNTVSTSHSGSFAPNPDLAAGAISLIVPGKSPTRKKS
 GPFGRSSAIGIENIQEVQEKRESPPAGQKTPDSGHVSQEPKSENSSTQSSPEMPTTKNRAETAAQRAE
 ALKDFSRSSSSASSFASVVEETEGVDGEDTGLSVSSSGTPHKRDSFIYSTWLEDSVSTTSGSSPGPSR
 SPHPDAGKLGDPACPEIKIQLEASEQHMPQLGC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6093_f04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_002885

ORF Size: 1989 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_002885.4</u>
RefSeq Size:	3334 bp
RefSeq ORF:	1992 bp
Locus ID:	5909
UniProt ID:	<u>P47736</u>
Cytogenetics:	1p36.12
MW:	73.3 kDa
Gene Summary:	This gene encodes a type of GTPase-activating-protein (GAP) that down-regulates the activity of the ras-related RAP1 protein. RAP1 acts as a molecular switch by cycling between an inactive GDP-bound form and an active GTP-bound form. The product of this gene, RAP1GAP, promotes the hydrolysis of bound GTP and hence returns RAP1 to the inactive state whereas other proteins, guanine nucleotide exchange factors (GEFs), act as RAP1 activators by facilitating the conversion of RAP1 from the GDP- to the GTP-bound form. In general, ras subfamily proteins, such as RAP1, play key roles in receptor-linked signaling pathways that control cell growth and differentiation. RAP1 plays a role in diverse processes such as cell proliferation, adhesion, differentiation, and embryogenesis. Alternative splicing results in multiple transcript variants encoding distinct proteins. [provided by RefSeq, Aug 2011]


Product images:



Circular map for RC208271



Western blot validation of overexpression lysate (Cat# [LY419032]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC208271 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified RAP1GAP protein (Cat# [TP308271]). The protein was produced from HEK293T cells transfected with RAP1GAP cDNA clone (Cat# RC208271) using MegaTran 2.0 (Cat# [TT210002]).