

Product datasheet for **RC208350**

KCNN3 (NM_170782) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KCNN3 (NM_170782) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	KCNN3
Synonyms:	hSK3; KCa2.3; SK3; SKCA3; ZLS3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC208350 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGAGAGACCTATAAAGGACTCCATGTTTTCGTTGGCCCTGAAATGCCTTATCAGTCTGTCCACCATCA
 TCCTTTTGGGCTTGATCATCGCCTACCACACACGTGAAGTCCAGCTCTTGTGATCGACAATGGCGCGGA
 TGACTGGCGGATAGCCATGACCTACGAGCGCATCCTGTACATCAGCCTGGAGATGCTGGTGTGCGCCATC
 CACCCATTCTGGCGAGTACAAGTTCTTCTGGACGGCACGCCTGGCCTTCTCTACACACCCTCCCGGG
 CGGAGGCCGATGTGGACATCATCTGTCTATCCCCATGTTCTGCGCCTGTACCTGATCGCCGAGTCAT
 GCTGCTGCACAGCAAGCTCTTACCAGTGCCTCGTCCCGCAGCATCGGGGCCCTCAACAAGATCACTTC
 AACACCCGCTTTGTCATGAAGACGCTCATGACCATCTGCCCTGGCACTGTGCTGCTGTTTCAGCATCT
 CTCTGTGGATCATTGCTGCTGGACCGTCCGTGTCTGTGAAAGGTACCATGACCAGCAGGACGTAAGTAG
 TAACTTTCTGGGTGCCATGTGGCTCATCTCCATCACATTCCTTTCCATTGGTTATGGGGACATGGTGCC
 CACACATACTGTGGGAAAGGTGTCTGTCTCCTCACTGGCATCATGGGTGCAAGGCTGCACTGCCCTTGTGG
 TGCCCGTGGTGGCCGAAAGCTGGAACACCAAAGCGGAGAAGCACGTTTATAACTTATGATGGACAC
 TCAGCTACCAAGCGGATCAAGAATGCTGCAGCCAATGTCCTTCGGGAAACATGGTTAATCTATAAACAC
 ACAAGCTGCTAAAGAAGATTGACCATGCCAAAGTGAGGAAACACCAGAGGAAGTTCTCCAAGCTATCC
 ACCAGTTGAGGAGCGTCAAGATGGAACAGAGGAAGCTGAGTGACCAAGCCAACACTCTGGTGGACCTTTC
 CAAGATGCAGAAATGTCATGTATGACTTAATCACAGAACTCAATGACCGGAGCGAAGACCTGGAGAAGCAG
 ATTGGCAGCCTGGAGTCAAGCTGGAGCATCTCACCGCCAGCTTCAACTCCCTGCCGTGCTCATCGCCG
 ACACCCGCGCCAGCAGCAGCAGCAGCTCTGTCTGCCATCATCGAGGCCCGGGGTGTGAGCGTGGCAGT
 GGGCACCAACCCACACCCCAATCTCCGATAGCCCCATTGGGGTCACTCCACCTCCTTCCCGACCCCGTAC
 ACAAGTTCAAGCAGTTGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC208350 protein sequence
 Red=Cloning site Green=Tags(s)

MERPIKDSMFLALKCLISLSTIILLGLIIAYHTREVQLFVIDNGADDWRIAMTYERILYISLEMLVCAI
 HPIPGEYKFFWTARLAFSYTPSRAEADVDIILSIPMFLRLYLARVMLLHSLKFTDASSRSIGALNKINF
 NTRFVMKTLMTICPGTVLLVFSISLWIIAAWTVRVCERYHDQDQVTSNFLGAMWLISITFLSIGYDMVP
 HTYCGKGVCLLTGIMGAGCTALVVAVVARKLELTKAEKHVHNFMMDTQLTKRIKNAANVLRWLIYKH
 TKLLKKIDHAKVRKHQRKFLQAIHQLRSVKMEQRKLSQANTLVDSLKMQNVMYDLITELNDRSEDLKQ
 IGSLESKLEHLTASFNSLPLLIADTLRQQQQLLSAIIIEARGVSVAVGTTHTPISDSPIGVSSTSFPPTY
 TSSSSC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6552_e09.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_170782

ORF Size: 1278 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_170782.3](#)

RefSeq Size: 11956 bp

RefSeq ORF: 1281 bp

Locus ID: 3782

UniProt ID: [Q9UGI6](#)


Cytogenetics: 1q21.3

Protein Families: Druggable Genome, Ion Channels: Potassium, Transmembrane


MW: 48.1 kDa

Gene Summary: Action potentials in vertebrate neurons are followed by an afterhyperpolarization (AHP) that may persist for several seconds and may have profound consequences for the firing pattern of the neuron. Each component of the AHP is kinetically distinct and is mediated by different calcium-activated potassium channels. This gene belongs to the KCNN family of potassium channels. It encodes an integral membrane protein that forms a voltage-independent calcium-activated channel, which is thought to regulate neuronal excitability by contributing to the slow component of synaptic AHP. This gene contains two CAG repeat regions in the coding sequence. It was thought that expansion of one or both of these repeats could lead to an increased susceptibility to schizophrenia or bipolar disorder, but studies indicate that this is probably not the case. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2011]

Product images:



Circular map for RC208350



Western blot validation of overexpression lysate (Cat# [LY406879]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC208350 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).