

## Product datasheet for **RC208355L4V**

### TRAPPC5 (NM\_174894) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TRAPPC5 (NM_174894) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRAPPC5
Synonyms:	TRS31
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_174894
ORF Size:	538 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208355).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_174894.1</a>
RefSeq Size:	814 bp
RefSeq ORF:	567 bp
Locus ID:	126003
UniProt ID:	<a href="#">Q8IURO</a>
Cytogenetics:	19p13.2
MW:	20.6 kDa



[View online »](#)

**Gene Summary:**

May play a role in vesicular transport from endoplasmic reticulum to Golgi.[UniProtKB/Swiss-Prot Function]