

Product datasheet for **RC208391A1V**

Human SFRP5 (NM_003015) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SFRP5 (NM_003015) AAV Particle
Tag:	Myc-DDK
Symbol:	SFRP5
Synonyms:	SARP3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208391 representing NM_003015 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCGGGCGGGCGGGCGGGGGGGCGTGCGGACGGCCGCGCTGGCGCTGCTGCTGGGGCGCTGCACT
GGGCGCCGGCGCGCTGCGAGGAGTACGACTACTATGGCTGGCAGGCCGAGCCGCTGCACGGCCGCTCCTA
CTCCAAGCCGCGCAGTGCCTTGACATCCCTGCCGACCTGCCGCTCTGCCACACGGTGGGCTACAAGCGC
ATGCGGCTGCCAACCCTGCTGGAGCACGAGAGCCTGGCCGAAGTGAAGCAGCAGGCGAGCAGCTGGCTGC
CGCTGCTGGCCAAGCGCTGCCACTCGGATACGACAGTCTTCTGTGCTCGCTCTTTCGCCCCGTCTGTCT
CGACCGGCCCATCTACCCGTGCCGCTCGCTGTGCGAGGCCGTGCGCGCCGGCTGCGCGCCGCTCATGGAG
GCCTACGGCTTCCCCTGGCCTGAGATGCTGCACTGCCACAAGTCCCCCTGGACAACGACCTCTGCATCG
CCGTGCAGTTCGGGCACCTGCCCGCCACCGCGCCTCCAGTGACCAAGATCTGCGCCAGTGTGAGATGGA
GCACAGTGCTGACGGCCTCATGGAGCAGATGTGCTCCAGTGACTTTGTGGTCAAATGCGCATCAAGGAG
ATCAAGATAGAGAATGGGGACCGGAAGCTGATTGGAGCCAGAAAAAGAAGCTGCTCAAGCCGGGCC
CCCTGAAGCGCAAGGACCAAGCGGCTGGTGTGCACATGAAGAATGGCGGGGCTGCCCTGCCACA
GCTGGACAGCCTGGCGGGCAGCTTCTGGTGTGCGCCGCAAGTGGATGGACAGCTGCTGCTCATGGCC
GTCTACCGCTGGGACAAGAAGAATAAGGAGATGAAGTTTGCAGTCAAATTCATGTTCTCCTACCCCTGCT
CCCTCTACTACCCTTTCTCTACGGGGCGGCAGAGCCAC

ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208391 representing NM_003015 Red=Cloning site Green=Tags(s) MRAAAAAGGVRTAALALLL GALHWAPARCEEYDYYGWQAEPLHGRSYSKPPQCLDIPADLPLCHTVGYKR MRLPNLLEHESLAEVKQQAASSWLPLLAKRCHSDTQVFLCSLFAPVCLDRPIYPCRS LCEAVRAGCAPLME AYGFPWPPEMLHCHKFPLDNDLCIAVQFGHLPATAPPVTKICAQCEMEHSADGLMEQMCS SDFVVKMRIKE IKIENGDRKLI GAQKKKKLLKPGPLKRKDTKRLVLHMKNGAGCPCPQLDSL AGSFLVMGRKVDGQ LLLMA VYRWDKKNKEMKFAVKFMFSYPCSLYYPFFYGAAEPH TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003015
ORF Size:	951 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003015.2 , NP_003006.1
RefSeq Size:	1905 bp
RefSeq ORF:	954 bp
Locus ID:	6425
UniProt ID:	Q5T4F7
Cytogenetics:	10q24.2
MW:	35.4 kDa