

Product datasheet for **RC208395L3V**

HOXC4 (NM_014620) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | HOXC4 (NM_014620) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | HOXC4 |
| Synonyms: | cp19; HOX3; HOX3E |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_014620 |
| ORF Size: | 792 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC208395). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_014620.4 |
| RefSeq Size: | 2328 bp |
| RefSeq ORF: | 795 bp |
| Locus ID: | 3221 |
| UniProt ID: | P09017 |
| Cytogenetics: | 12q13.13 |
| Protein Families: | Transcription Factors |
| MW: | 29.8 kDa |



[View online »](#)

Gene Summary:

This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene, HOXC4, is one of several homeobox HOXC genes located in a cluster on chromosome 12. Three genes, HOXC5, HOXC4 and HOXC6, share a 5' non-coding exon. Transcripts may include the shared exon spliced to the gene-specific exons, or they may include only the gene-specific exons. Two alternatively spliced variants that encode the same protein have been described for HOXC4. Transcript variant one includes the shared exon, and transcript variant two includes only gene-specific exons. [provided by RefSeq, Jul 2008]