

Product datasheet for **RC208398L4V**

Munc18c (STXBP3) (NM_007269) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Munc18c (STXBP3) (NM_007269) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Munc18c |
| Synonyms: | MUNC18-3; MUNC18C; PSP; UNC-18C |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_007269 |
| ORF Size: | 1776 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC208398). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_007269.2 |
| RefSeq Size: | 2522 bp |
| RefSeq ORF: | 1779 bp |
| Locus ID: | 6814 |
| UniProt ID: | O00186 |
| Cytogenetics: | 1p13.3 |
| Domains: | Sec1 |
| Protein Families: | Stem cell - Pluripotency |



[View online »](#)

MW: 67.8 kDa

Gene Summary: Together with STX4 and VAMP2, may play a role in insulin-dependent movement of GLUT4 and in docking/fusion of intracellular GLUT4-containing vesicles with the cell surface in adipocytes.[UniProtKB/Swiss-Prot Function]