

Product datasheet for **RC208410A1V**

Human EMP1 (NM_001423) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human EMP1 (NM_001423) AAV Particle
Tag: Myc-DDK
Symbol: EMP1
Synonyms: CL-20; EMP-1; TMP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208410 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTTGGTATTGCTGGCTGGTATCTTTGTGGTCCACATCGCTACTGTTATTATGCTATTTGTTAGCACCA
TTGCCAATGTCTGGTTGGTTTCCAATACGGTAGATGCATCAGTAGGTCTTTGGAAAACTGTACCAACAT
TAGCTGCAGTGACAGCCTGTCTATGCCAGTGAAGATGCCCTCAAGACAGTGCAGGCCTTCATGATTCTC
TCTATCATCTTCTGTGTCTATGCCCTCCTGGTCTTCGTGTTCCAGCTCTTCACCATGGAGAAGGGAAACC
GGTTCCTCCTCAGGGGCCACACACTGGTGTGCTGGCTGTGCATTCTTGTGGGGGTGCCATCTACAC
TAGTCATTATGCGAATCGTGATGGAACGAGTATCACCACGGCTATTCTACATCCTGGGCTGGATCTGC
TTCTGCTTCAGCTTCATCATCGGCGTTCTCTATCTGGTCTGAGAAAGAAA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208410 protein sequence
Red=Cloning site Green=Tags(s)

MLVLLAGIFVVIATVIMLFVSTIANVWLVSNTVDASVGLWKNCTNISCSDSLSEYSEDALKTVQAFMIL
SIIFCVIALLVFVQLFTMEKGNRFFLSGATTLVCWLCILVGVSITYSHYANRDGTQYHHGYSYILGWIC
FCFSFIIGVLYLVLRKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_001423
ORF Size:	471 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001423.1</u>
RefSeq Size:	2804 bp
RefSeq ORF:	474 bp
Locus ID:	2012
UniProt ID:	<u>P54849</u>
Cytogenetics:	12p13.1
MW:	17.6 kDa