

Product datasheet for **RC208433A1V**

Human PLEKHA3 (NM_019091) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PLEKHA3 (NM_019091) AAV Particle
Tag:	Myc-DDK
Symbol:	PLEKHA3
Synonyms:	FAPP1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208433 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGGGGGTGTGTACAAGTGGACCAACTATCTCACAGGCTGGCAGCCTCGTTGGTTTGTTTAGATA
ATGGAATCTTATCCTACTATGATTCACAAGATGATGTTTGCAAAGGGAGCAAAGGAAGCATAAAGATGGC
AGTTTGTGAAATTAAGTTCATTCAGCAGACAACACAAGAATGGAATTAATCATTCTGGAGAGCAGCAT
TTCTACATGAAGGCAGTGAATGCAGCTGAAAGACAGAGGTGGCTGGTCTGCTCTGGGGAGCTCCAAAGCAT
GTTTGACTGATACAAGGACTAAAAAGAAAAAGAAATAAGTGAAACCAAGTGAATCGCTGAAAACCAAAAT
GTCTGAACCTTCGCCTCTACTGTGACCTCTTAATGCAGCAAGTTCATACAATACAGGAATTTGTTCAACCAT
GATGAGAATCATTTCATCTCCTAGTGCAGAGAACATGAATGAAGCCTCTTCTCTGCTTAGTGCCACGTGTA
ACACATTCATCACAACGCTTGAGGAATGTGTGAAGATAGCCAATGCCAAGTTTAAACCTGAGATGTTTCA
ACTGCACCATCCGGATCCCTTAGTTTCTCCTGTGTACCTTCTCCTGTTCAAATGATGAAGCGTTCTGTCT
AGCCACCCTGGTCTTGCAGTTCAGAGAGGAGTAGCCACTCTATAAAAGAACCAGTATCTACACTTCACC
GACTCTCCAGCGACGCCGAAGAACCTACTCAGATACAGATTCTTGATGATATTCTCTTGAAGACCC
AGATAGACCTGTTCACTGTTCAAAAAATACACTTAATGGAGATTGGCATCAGCAACCATTCTGAAGAA
AGCAGACTTATGGCCAAAAACAATCTGAATCAGAAGATACTCTCCATCCTTCTCTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208433 protein sequence Red=Cloning site Green=Tags(s) MEGVLYKWTNYLTGWQPRWFLDNGILSYYSQDDVCKGSGSIKMAVCEIKVHSADNTRMELIIPGEQH FYMKAVNAAERQRLVALGSSKACLTDTRTKKEKEISETSESLKTKMSELRLYCDLLMQQVHTIQEFVHH DENHSSPSAENMNEASSLLSATCNTFITLLEECVKIANAKFKPEMFQLHHPDPLVSPVSPSPVQMMKRSV SHPGSCSSERSHSIPEVSTLHRLSQRRRTYSDTDSCSDIPLDPDRPVHCSKNTLNGDLASATIPPEE SRLMAKKQSEEDTLPSFSS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_019091
ORF Size:	900 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_019091.3</u>
RefSeq Size:	2516 bp
RefSeq ORF:	903 bp
Locus ID:	65977
UniProt ID:	<u>Q9HB20</u>
Cytogenetics:	2q31.2
MW:	33.9 kDa