

## Product datasheet for **RC208438L2V**

### SMARCAD1 (NM\_020159) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SMARCAD1 (NM_020159) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SMARCAD1
Synonyms:	ADERM; BASNS; ETL1; HEL1; HRZ
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_020159
ORF Size:	3084 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208438).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020159.1</a>
RefSeq Size:	5128 bp
RefSeq ORF:	3081 bp
Locus ID:	56916
UniProt ID:	<a href="#">Q9H4L7</a>
Cytogenetics:	4q22.3
Domains:	SNF2_N, DEAD, helicase_C, CUE
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Transcription Factors



[View online »](#)

**MW:** 117.6 kDa

**Gene Summary:** This gene encodes a member of the SNF subfamily of helicase proteins. The encoded protein plays a critical role in the restoration of heterochromatin organization and propagation of epigenetic patterns following DNA replication by mediating histone H3/H4 deacetylation. Mutations in this gene are associated with adermatoglyphia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]