

Product datasheet for **RC208448L3V**

WARS2 (NM_015836) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	WARS2 (NM_015836) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WARS2
Synonyms:	mtTrpRS; NEMMLAS; TrpRS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015836
ORF Size:	1080 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208448).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015836.3
RefSeq Size:	2806 bp
RefSeq ORF:	1083 bp
Locus ID:	10352
UniProt ID:	Q9UGM6
Cytogenetics:	1p12
Domains:	tRNA-synt_1b
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Aminoacyl-tRNA biosynthesis, Tryptophan metabolism

MW: 40.2 kDa

Gene Summary: Aminoacyl-tRNA synthetases catalyze the aminoacylation of tRNA by their cognate amino acid. Because of their central role in linking amino acids with nucleotide triplets contained in tRNAs, aminoacyl-tRNA synthetases are thought to be among the first proteins that appeared in evolution. Two forms of tryptophanyl-tRNA synthetase exist, a cytoplasmic form, named WARS, and a mitochondrial form, named WARS2. This gene encodes the mitochondrial tryptophanyl-tRNA synthetase. Two alternative transcripts encoding different isoforms have been described. [provided by RefSeq, Jul 2008]