

Product datasheet for **RC208458A1V**

Human FBXO4 (NM_012176) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FBXO4 (NM_012176) AAV Particle
Tag:	Myc-DDK
Symbol:	FBXO4
Synonyms:	FBX4
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208458 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGGAAGCGAGCCGCGCAGCGGAACAACTCGCCGCCGCCCTTCAGCGACTGGGGCCGCCTGG
AGGCGGCCATCCTCAGCGGCTGGAAGACCTTCTGGCAGTCAGTGAGCAAGGAGAGGGTGGCGCTACGAC
CTCAGGGGAGGAGGTGGATGAGCGGCCAGCACCTGACGCGGCTGCCGATTGATGTACAGCTATATATT
TTGTCCTTTCTTCACCTCATGATCTGTGTCAGTTGGGAAGTACAAATCATTATTGGAATGAACTGTAA
GAGATCCAATTCTGTGGAGATACTTTTTGTTGAGGGATCTTCTTCTTGGTCTTCTGTTGACTGGAAGTC
TCTTCCAGATCTAGAAATCTTAAAAAGCCTATATCTGAGGTCACTGATGGTGCATTTTTTGACTACATG
GCAGTCTATAGAATGTGCTGTCCATACACAAGAAGAGCTTCAAAATCCAGCCGTCTATGTATGGAGCTG
TCACTTCTTTTTTACACTCCCTGATCATTGAGATGAACACGATTTGCTATGTTTGGACCAGGTTTGA
AGAATTGAATACCTCTTTGGTGTGAGCTTGATGTCTCAGAGGAACCTTGCCCAACAGCTGGTTTGCCT
CAGAGGCAGATTGATGGTATTGGATCAGGAGTCAATTTTCAGTTGAACAACCAACATAAATTCAACATTC
TAATCTTATATTCACTACCAGAAAGGAAAGAGATAGAGCAAGGGAAGAGCATAAAGTGCAGTTAACA
GATGTTCACTGACACAATGAAGGTGATGATCAACAAGGAAGCCGGTACAGTGTGATTCCACAGATTCAA
AAAGTGTGTGAAGTTGTAGATGGTTCATCTATGTTGCAAATGCTGAAGCTCATAAAAGACATGAATGCC
AAGATGAATTTTCTCATATTATGGCAATGACAGATCCAGCCTTTGGGTCTTCGGGAAGACCATTGTTGGT
TTTATCTTGTATTTCTCAAGGGGATGTAAAAAGATGCCTGTTTTTATTGGCTCATGAGCTGCATCTG
AATCTTCTAAATCACCATGGCTGGTCCAGGATACAGAGGCTGAAACTCTGACTGGTTTTTTGAATGGCA
TTGAGTGGATTCTTGAAGAAGTGAATCTAAGCGTGCAAGA

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208458 protein sequence Red=Cloning site Green=Tags(s) MAGSEPRSGTNSPPPPFSDWGRLEAAILSGWKTFWQSVSKERVARTTSREEVDEAASTLTRLPIDVQLYI LSFLSPHDL CQLGSTNHYYNETVRDPI LWRYFLLRDLPSWSSVDWKS L PDL EILKKPISEVTDGAFFDYM AVYRMCCPYTRASKSSRP MYGAVTSFLHSLIIQNEPRFAMFGPGL EELNTSLVLSLMSSEELCPTAGLP QRQIDGIGSGVNFQLNNQHKFNILILYSTTRKERDRAREEHTSAVNKMFSRHNEGDDQQGSRYSVIPQIQ KVCEVVDGFIYVANA EAHKRHEWQDEF SHIMAMTDPAFGSSGRPLLVLSCISQGDVYKRMPCFYLAHELHL NLLNHPWL VQDTEAETLTGFLNGIEWILEEVESKRAR SGP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_012176
ORF Size:	1161 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_012176.2
RefSeq Size:	1523 bp
RefSeq ORF:	1164 bp
Locus ID:	26272
UniProt ID:	Q9UKT5
Cytogenetics:	5p13.1
MW:	44.1 kDa