

Product datasheet for **RC208465L3V**

CEL3 (NM_007185) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CEL3 (NM_007185) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CEL3
Synonyms:	BRUNOL1; CAGH4; ERDA4; ETR-1; TNRC4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007185
ORF Size:	1035 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208465).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007185.3
RefSeq Size:	5615 bp
RefSeq ORF:	1398 bp
Locus ID:	11189
UniProt ID:	Q5SZQ8
Cytogenetics:	1q21.3
Domains:	RRM
Protein Families:	Transcription Factors



[View online »](#)

MW: 50.5 kDa

Gene Summary: Members of the CELFB/BRUNOL protein family contain two N-terminal RNA recognition motif (RRM) domains, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RRM domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involved in mRNA editing, and translation. Multiple alternatively spliced transcript variants encoding different isoforms have been identified in this gene. [provided by RefSeq, Feb 2010]