

Product datasheet for **RC208481L4V**

APOL6 (NM_030641) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	APOL6 (NM_030641) Human Tagged ORF Clone Lentiviral Particle
Symbol:	APOL6
Synonyms:	APOL-VI; APOLVI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030641
ORF Size:	1029 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208481).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030641.3
RefSeq Size:	10156 bp
RefSeq ORF:	1032 bp
Locus ID:	80830
UniProt ID:	Q9BWW8
Cytogenetics:	22q12.3
MW:	38.1 kDa



[View online »](#)

Gene Summary:

This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to organelles. [provided by RefSeq, Jul 2008]