

Product datasheet for **RC208511L1V**

FAM8A1 (NM_016255) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM8A1 (NM_016255) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM8A1
Synonyms:	AHCP
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016255
ORF Size:	1239 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208511).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016255.2
RefSeq Size:	4762 bp
RefSeq ORF:	1242 bp
Locus ID:	51439
UniProt ID:	Q9UBU6
Cytogenetics:	6p22.3
Protein Families:	Transmembrane
MW:	44.1 kDa



[View online »](#)

Gene Summary:

Plays a role in the assembly of the HRD1 complex, a complex involved in the ubiquitin-proteasome-dependent process of ER-associated degradation (ERAD).[UniProtKB/Swiss-Prot Function]