

Product datasheet for **RC208515A1V**

Human PILRB (NM_013440) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PILRB (NM_013440) AAV Particle
Tag: Myc-DDK
Symbol: PILRB
Synonyms: FDFACT1; FDFACT2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208515 ORF sequence
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGTGCGCCCTGCTGCTGCCCTGCTGCTCCTGCTGCAGCCGCCAGCATTTCTGCAGCCTGGTGGCT
CCACAGGATCTGGTCCAAGCTACCTTTATGGGGTCACTCAACCAAAACACCTCTCAGCCTCCATGGGTG
CTCTGTGGAATCCCTTCTCCTTCTATTACCCCTGGGAGTTAGCCATAGTTCCCAACGTGAGAATATCC
TGGAGACGGGGCACTTCCACGGGCAGTCCTTCTACAGCACAAGGCCGCTTCCATTACAAGGATTATG
TGAACCGGCTCTTTCTGAAGTGGACAGAGGGTCAGGAGAGCGGCTTCTCAGGATCTCAAACCTGCGGAA
GGAGGACCAAGTCTGTGTATTTCTGCCGAGTCGAGCTGGACACCCGGAGATCAGGGAGGCAGCAGTTGCAG
TCCATCAAGGGGACCAAACTCACCATCACCCAGGCTGTCAACAACCACCACCTGGAGGCCAGCAGCA
CAACCACCATAGCCGGCCTCAGGGTCACAGAAAGCAAAGGGCACTCAGAATCATGGCACCTAAGTCTGGA
CACTGCCATCAGGGTTGATTGGCTGTCGCTGTGCTCAAACTGTCAATTTGGGACTGCTGTGCCTCCTC
CTCCTGTGGTGGAGGAGAAGAAAGGTAGCAGGGCGCCAAGCAGTGACTTC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208515 protein sequence
Red=Cloning site **Green**=Tags(s)

MGRPLLLPLLLLQPPAFLQPGSGTSGSPSYLYGVTQPKHLSASMGGSVETPFYFPWELAIIVPNVRIS
WRRGHFHGQSFYSTRPPSIHKDYVNRFLNWTEGQESGFLRISNLRKEDQSVYFCRVELDTRRSGRQQLQ
SIKGTKLITQAVTTTTTWRPSSTTTIAGLRVTESKGHSESWHLSLDTAIRVALAVAVLKTIVILGLLCLL
LLWRRRKGSRAPSSDF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_013440
ORF Size:	681 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_013440.3</u> , <u>NP_038468.3</u>
RefSeq Size:	3632 bp
RefSeq ORF:	683 bp
Locus ID:	29990
Cytogenetics:	7q22.1
MW:	25.5 kDa