

Product datasheet for **RC208527L4V**

TIMELESS (NM_003920) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TIMELESS (NM_003920) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TIMELESS
Synonyms:	hTIM; TIM; TIM1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003920
ORF Size:	3624 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208527).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003920.2
RefSeq Size:	5154 bp
RefSeq ORF:	3627 bp
Locus ID:	8914
UniProt ID:	Q9UNS1
Cytogenetics:	12q13.3
Domains:	TIMELESS, TIMELESS_C
Protein Families:	Transcription Factors



[View online »](#)

MW: 138.7 kDa

Gene Summary: The protein encoded by this gene is highly conserved and is involved in cell survival after damage or stress, increase in DNA polymerase epsilon activity, maintenance of telomere length, and epithelial cell morphogenesis. The encoded protein also plays a role in the circadian rhythm autoregulatory loop, interacting with the PERIOD genes (PER1, PER2, and PER3) and others to downregulate activation of PER1 by CLOCK/ARNTL. Changes in this gene or its expression may promote prostate cancer, lung cancer, breast cancer, and mental disorders. [provided by RefSeq, Feb 2014]