

Product datasheet for **RC208536L4V**

ADSSL1 (ADSS1) (NM_152328) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ADSSL1 (ADSS1) (NM_152328) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADSS1
Synonyms:	ADSSL1; MPD5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152328
ORF Size:	1371 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208536).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152328.3
RefSeq Size:	1769 bp
RefSeq ORF:	1374 bp
Locus ID:	122622
UniProt ID:	Q8N142
Cytogenetics:	14q32.33
Domains:	Adenylsucc_synt
Protein Pathways:	Alanine, aspartate and glutamate metabolism, Metabolic pathways, Purine metabolism



[View online »](#)

MW: 50.2 kDa

Gene Summary: This gene encodes a member of the adenylosuccinate synthase family of proteins. The encoded muscle-specific enzyme plays a role in the purine nucleotide cycle by catalyzing the first step in the conversion of inosine monophosphate (IMP) to adenosine monophosphate (AMP). Mutations in this gene may cause adolescent onset distal myopathy. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]