

## Product datasheet for **RC208546L4V**

### ZNF179 (RNF112) (NM\_007148) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF179 (RNF112) (NM_007148) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF179
Synonyms:	BFP; ZNF179
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007148
ORF Size:	1893 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208546).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007148.2</a>
RefSeq Size:	3101 bp
RefSeq ORF:	1896 bp
Locus ID:	7732
UniProt ID:	<a href="#">Q9ULX5</a>
Cytogenetics:	17p11.2
Protein Families:	Transmembrane
MW:	68.1 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a member of the RING finger protein family of transcription factors. The protein is primarily expressed in brain. The gene is located within the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008]