

Product datasheet for **RC208547A1V**

Human DcR2 (TNFRSF10D) (NM_003840) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DcR2 (TNFRSF10D) (NM_003840) AAV Particle
Tag:	Myc-DDK
Symbol:	DcR2
Synonyms:	CD264; DCR2; TRAIL-R4; TRAILR4; TRUND
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208547 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGGACTTTGGGGACAAAGCGTCCCGACCGCTCGAGCGCTCGAGCAGGGCGCTATCCAGGAGCCAGGA
CAGCGTCGGGAACAGACCATGGCTCCTGGACCCCAAGATCCTTAAGTTCGTCGCTTTCATCGTCGCGGT
TCTGCTGCCGGTCCGGTTGACTCTGCCACCATCCCCCGGCAGGACGAAGTTCAGCAGACAGTGGCC
CCACAGCAACAGAGGCGCAGCCTCAAGGAGGAGGAGTGTCAGCAGGATCTCATAGATCAGAATATACTG
GAGCCTGTAACCCGTGCACAGAGGGTGTGGATTACACCATGCTTCCAACAATTTGCCTTCTTGCTGCT
ATGTACAGTTTGTAAATCAGGTCAAACAATAAAAGTTCTGTACACGACGACAGAGACACCGTGTGTGAG
TGTGAAAAAGGAAGCTTCCAGGATAAAAACTCCCTGAGATGTGCCGGACGTGTAGAACAGGGTGTCCCA
GAGGGATGGTCAAGGTCAAGTCAATTGTACGCCCCGAGTGACATCAAGTGCAAAATGAATCAGCTGCCAG
TTCCACTGGGAAAACCCAGCAGCGGAGGAGACAGTGACCACCATCCTGGGGATGCTTGCTCTCCCTAT
CACTACCTTATCATCATAGTGGTTTTAGTCATCATTTTAGCTGTGGTTGTGGTTGGCTTTTCATGTCGGA
AGAAATTCATTTCTTACCTCAAAGGCATCTGCTCAGGTGGTGGAGGAGTCCGAACGTGTGCACAGAGT
CCTTTTCCGGCGCGTTCATGTCCTTCACGAGTTCCTGGGCGGAGGACAATGCCCGCAACGAGACCTG
AGTAACAGATACTTGACGCCACCCAGGTCTCTGAGCAGGAAATCCAAGTTCAGGAGCTGGCAGAGCTAA
CAGGTGTGACTGTAGAGTCGCCAGAGGAGCCACAGCGTCTGCTGGAACAGGCAGAAGCTGAAGGGTGTCA
GAGGAGGAGGCTGCTGGTTCCAGTGAATGACGCTGACTCCGCTGACATCAGCACCTTGCTGGATGCCTCG
GCAACACTGGAAGAAGGACATGCAAAGGAAACAATTCAGGACCAACTGGTGGGCTCCGAAAAGCTCTTTT
ATGAAGAAGATGAGGCAGGCTCTGCTACGTCCTGCCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208547 protein sequence Red=Cloning site Green=Tags(s) MGLWGQSVPTASSARAGRYPGARTASGTRPWLLDPKILKFVVFIVAVLLPVRVDSATIPRQDEVPPQQTVA PQQRRSLKEECPAGSHRSEYTGACNPCTEGVDYTIASNNLPSCLLCTVCKSGQTNKSSCTTTTRDVCQ CEKGSFQDKNSPEMCRCTGTCPRGMVKVSNCTPRSDIKCKNESAASSTGKTPAAEETVTTILGMLASPY HYLIIIVVLVIIILAVVVVGFSCRKKFISYLKGICSGGGGGPERVHRVLFRRRSCPSRVPGAEDNARNETL SNRYLQPTQVSEQEIQQQELAEITGVTVESPEEPQRLLEQAEAGCQRRRLLPVNDADSADISTLLDAS ATLEEGHAKETIQDQLVGSEKLFYEDEAGSATSC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003840
ORF Size:	1158 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003840.3
RefSeq Size:	3569 bp
RefSeq ORF:	1161 bp
Locus ID:	8793
UniProt ID:	Q9UBN6
Cytogenetics:	8p21.3
MW:	41.8 kDa