

Product datasheet for **RC208548L1V**

PPP1R16A (NM_032902) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPP1R16A (NM_032902) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PPP1R16A
Synonyms:	MYPT3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_032902
ORF Size:	1584 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208548).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032902.5
RefSeq Size:	2326 bp
RefSeq ORF:	1587 bp
Locus ID:	84988
UniProt ID:	Q96I34
Cytogenetics:	8q24.3
Domains:	ANK
MW:	57.6 kDa



[View online »](#)

Gene Summary:

Myosin light chain kinase and phosphatase (MLCP) complexes control the phosphorylation states of regulatory myosin light chains, which is crucial for muscle and intracellular movement. MLCPs typically contain a catalytic protein phosphatase 1 (PP1c) subunit, a myosin phosphatase targeting (MYPT) subunit, and another smaller subunit. The protein encoded by this gene represents an MYPT subunit, which is responsible for directing PP1c to its intended targets. However, while the phosphorylation of other MYPT members results in PP1c inactivation, phosphorylation of the encoded protein by protein kinase A results in PP1c activation. [provided by RefSeq, Jan 2020]