

## Product datasheet for **RC208556L2V**

### Heparanase 1 (HPSE) (NM\_006665) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Heparanase 1 (HPSE) (NM_006665) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Heparanase 1
Synonyms:	HPA; HPA1; HPR1; HPSE1; HSE1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006665
ORF Size:	1629 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208556).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006665.3</a> , <a href="#">NP_006656.1</a>
RefSeq Size:	4668 bp
RefSeq ORF:	1632 bp
Locus ID:	10855
UniProt ID:	<a href="#">Q9Y251</a>
Cytogenetics:	4q21.23
Domains:	Glyco_hydro_79n
Protein Families:	Secreted Protein



[View online »](#)

**Protein Pathways:** Glycosaminoglycan degradation, Metabolic pathways

**MW:** 61.2 kDa

**Gene Summary:** Heparan sulfate proteoglycans are major components of the basement membrane and extracellular matrix. The protein encoded by this gene is an enzyme that cleaves heparan sulfate proteoglycans to permit cell movement through remodeling of the extracellular matrix. In addition, this cleavage can release bioactive molecules from the extracellular matrix. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]