

Product datasheet for RC208557L4

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

CD163 (NM_004244) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: CD163 (NM_004244) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: CD163

Synonyms: M130; MM130; SCARI1

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)


E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC208557).

Sequence:

Restriction Sites: Sgfl-Rsrll

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_004244

ORF Size: 3468 bp





CD163 (NM_004244) Human Tagged Lenti ORF Clone - RC208557L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 004244.3</u>

RefSeq Size:4190 bpRefSeq ORF:3471 bpLocus ID:9332

UniProt ID: Q86VB7

Cytogenetics: 12p13.31

Domains: SR

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

MW: 125.4 kDa

Gene Summary: The protein encoded by this gene is a member of the scavenger receptor cysteine-rich (SRCR)

superfamily, and is exclusively expressed in monocytes and macrophages. It functions as an


acute phase-regulated receptor involved in the clearance and endocytosis of

hemoglobin/haptoglobin complexes by macrophages, and may thereby protect tissues from free hemoglobin-mediated oxidative damage. This protein may also function as an innate immune sensor for bacteria and inducer of local inflammation. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq,


Aug 2011]



Product images:



Circular map for RC208557L4



Double digestion of RC208557L4 using Sgfl and Rsrll