

Product datasheet for **RC208560A1V**

Human CDK6 (NM_001259) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CDK6 (NM_001259) AAV Particle
Tag:	Myc-DDK
Symbol:	CDK6
Synonyms:	MCPH12; PLSTIRE
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208560 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGAAGGACGGCCTGTGCCGCGCTGACCAGCAGTACGAATGCGTGGCGGAGATCGGGGAGGGCGCCT
 ATGGGAAGGTGTTCAAGGCCGCGCACTGAAGAACGGAGGCCGTTTCGTGGCGTTGAAGCGCGTGCGGGT
 GCAGACCGCGGAGGAGGGCATGCCGCTCTCCACCATCCGCGAGGTGGCGGTGCTGAGGCACCTGGAGACC
 TTCGAGCACCCCAACGTGGTCAGGTTGTTTGATGTGTGCACAGTGTACGAACAGACAGAGAAACCAAC
 TAACCTTTAGTGTGTTGAACATGTCGATCAAGACTTGACCCTTACTTGGATAAAGTTCAGAGCCTGGAGT
 GCCCACTGAAACCATAAAGGATATGATGTTTCAGCTTCTCCGAGGTCTGGACTTTCTTCATTACACCGA
 GTAGTGCATCGCGATCTAAAACACAGAACATTCTGGTGACCAGCAGCGGACAAATAAACTCGTGACT
 TCGGCCTTGCCCGCATCTATAGTTTCCAGATGGCTCTAACCTCAGTGGTCGTACGCTGTGGTACAGAGC
 ACCCGAAGTCTTGCTCCAGTCCAGCTACGCCACCCCGTGGATCTCTGGAGTGTGGCTGCATATTTGCA
 GAAATGTTTCGTAGAAAGCCTCTTTTCGTGGAAGTTCAGATGTTGATCAACTAGGAAAAATCTTGACG
 TGATTGGACTCCAGGAGAAGAAGACTGGCCTAGAGATGTTGCCCTTCCAGGCAGGCTTTTCATTCAA
 ATCTGCCCAACCAATTGAGAAGTTGTAAACAGATATCGATGAAGTACGAAAGACCTACTTCTGAAGTGT
 TTGACATTTAACCAGCCAAAAGAATATCTGCCTACAGTGCCCTGTCTACCCATACTTCCAGGACCTGG
 AAAGGTGCAAAGAAAACCTGGATTCCACCTGCCGCCAGCCAGAACACCTCGGAGCTGAATACAGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC208560 protein sequence Red=Cloning site Green=Tags(s) MEKDGLCRADQQYECVAEIGEGAYGKVFKARDLKNNGRFVALKRVRVQTGEEGMPLSTIREVAVLRHLET FEHPNVVRLFDVCTVSRTDRETKLTLVFEHVDQDLTTYLDKVPEPGVPTETIKDMMFQLLRGLDFLHSHR VVHRDLKPQNILVTSSGQIKLADFGGLARIYSFQMALTSVVVTLWYRAPEVLLQSSYATPVDLWSVGCIFA EMFRRKPLFRGSSDQDLGKILDVIGLPGEEDWPRDVALPRQAFHKSQAQPIEFVTDIDELGKDLLLK LTFNPAKRISAYSALSHPYFQDLERCKENLDShLPPSQNTSELNTA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001259
ORF Size:	978 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001259.5</u>
RefSeq Size:	11628 bp
RefSeq ORF:	981 bp
Locus ID:	1021
UniProt ID:	<u>Q00534</u>
Cytogenetics:	7q21.2
MW:	36.9 kDa