

## Product datasheet for **RC208570L1V**

### SLC2A4RG (NM\_020062) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC2A4RG (NM_020062) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC2A4RG
Synonyms:	GEF; HDBP-1; HDBP1; Si-1-2; Si-1-2-19
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_020062
ORF Size:	1161 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208570).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020062.3</a>
RefSeq Size:	2296 bp
RefSeq ORF:	1164 bp
Locus ID:	56731
UniProt ID:	<a href="#">Q9NR83</a>
Cytogenetics:	20q13.33
Domains:	zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 41.1 kDa

**Gene Summary:** The protein encoded by this gene is a nuclear transcription factor involved in the activation of the solute carrier family 2 member 4 gene. The encoded protein interacts with another transcription factor, myocyte enhancer factor 2, to activate transcription of this gene. [provided by RefSeq, Jul 2008]