

Product datasheet for **RC208577**

VRTN (NM_018228) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VRTN (NM_018228) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	VRTN
Synonyms:	C14orf115; vertnin
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC208577 ORF sequence

Red=Cloning site Blue=ORF Green=Tags(s)

 TTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

 ATGACTTCTCGGAACCAGCTGGTGCAGAAGGTGCTGCAGGAGCTGCAGGAAGCAGTGAGGTGCGAAGGCC
TGGAGGGTCTCATAGGTGCTTCCTTGGAGGCCAAGCAGGTCTGTCTTCCTTCACTCTCCCCACCTGCCG
GGAGGGAGGCCCTGGCCTCCAGGTGCTGGAAGTGGACTCGGTGGCCCTGAGCCTGTATCCAGAAGATGCT
CCACGGAACATGCTGCCGCTGGTGTGCAAGGGGAGGGCAGCCTGCTGTTGAGGCGGCCAGCATGCTGC
TGTGGGGTGACGAGGCCCTCAGCCTGGAGCTGCGGGCCCGCACCGTGGTAGAGATGCTGTCACAGACA
CTACTACCTCCAGGCATGATCGACTCCAAAGTATGCTGCAGGCCGTGCGCTACTCCCTATGCTCTGAG
GAGTCCCTGAGATGACCAGCTTGCCCCCGCCACGCTGGAGGCCATCTTCGATGCCGACGTCAGGCCCT
CCTGTTTCCCGCAGCAGCTTCTCAACGTGTGGCACTTGTATGCTCTCGCCTCTGCTCTCCAGCGGAACAT
CTACTCCATCTACCCATGCGCAACCTCAAGATCCGGCCCTACTTCAACCGTGTATCCGGCCCCCGCCG
TGCAGACCACTGCGCTCCACGCTGCACATCATGTGGGCTGGCCAGCCCCCTACCAGCCACTTCTTCCGCC
ACCACTACTTTGCCCCCTGTGGTGGGGCTGGAAGAGGTGGAGGCTGAAGGTGCCCTTGGCGTGGCCCCAGC
TCTTCCAGCCCTGGCCCCACTCTCATCGCCGGCCAAGACCCTGGAGCTGCTCAACCGTGAACCTGGCCTC
AGCTACTCTCACCTCTGTGAGCGCTACAGCGTCACCAAAAGCACCTTCTACCGCTGGCGGGCAGTCCC
AGGAGCACCGGCAGAAGTTGCTGCCCGCTTCTCCGCCAAGCACTTCTGCAGGACAGCTTCCACCGGGG
GGGCGTCGTGCCACTTCAGCAGTTCTCCAGCGGTTCCCGGAGATCTCCCGCTCAACCTACTATGCTGTG
AAGCATGAGCTGCTGGGCTCTGGCACCTGCCCGGCTTGGCCCCAGGGAGGTGCTGGGCATGGAGGAGC
TAGAGAAGCTGCCGAGGAGCAGGTGGCTGAGGAGGAGCTGGAGTGTCCGCACTGGCGGTGTCAAGCCC
TGGAAATGGTCTTAATGCAGCGGGCCAAGTTGTACCTGGAGCATTGCATCTCCCTGAACACACTGGTACCC
TATCGCTGCTTCAAACGCAAGTTCCCTGGCATCTCACGGTCCACTTATTATAATTGGCGGCGAAAGGCC
TCCGGAGGAACCCAGCTTCAAGCCGGCACCAGCCCTCTCTGCTGCTGGGACTCCCGAGCTAGCATCTGT
TGGGGAAGGGGCTGTAATTCCTTGAAGAGTGAGGCGGAAGAGGGGGCAGGGAATGCCACAGGTGAGGAC
CCTCCCGCCCCGGGAGCTCTGCCACTAAGGATGCCCTGTCCCGTTGGCAGAGGCGTCTGCGCAGGG
CTGCCCCGAGGCAGGTGCTGAGTGGGCATCTCCCTTTCTGCCGCTTCCGCTCCGCTACCCAGCCTGTC
ACCTTCTGCCTTTTGGGTCTGGAAGAGTCTTGCTCGGGTTGGCCAGAGGCCTGTCAAACCTTCAGGTG
CCGGTCCCGACCTTGGGCAAAGGGGGCAGGAGGCTGAGGAGAAGCAGGAGAAGGAGGCTGGCAGGGATG
TGACAGCTGTGATGGCCCCACCTGTGGGGCTTCTTCAGAAGATGTAGAGGGAGGGCCTTCCAGAGAGGG
GGCCCTGCAGGAGGGGGCCACAGCCAGGGCCAGCCCCACAGTGGGCCCTTGTGAGCCAACCTGTGGTG
GCAGCAGCGGGTGGCAGGGATGGCCGATGCTGGTGATGGACATGATCGTACCACGAAGTTCAAGGCC
AGGCCAAGCTGTTCTTGAGAAGCGCTTCCAGTCCAAGAGCTTCCCTCTACAAGGAGTTCAGTGCCCT
CTTCCCTCACTGCCGCTCCACATACTACATGTGGAAGCGAGCCCTCTATGACGGCTGACCCTGGTA
GATGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208577 protein sequence
 Red=Cloning site Green=Tags(s)


MTSRNQLVQKVLQELQEAVECEGLEGLIGASLEAKQVLSSFTLPTCREGGPGLQVLEVDVSVALSLYPEDA
 PRNMLPLVCKGEGSLLFEAASMLLWGDAGLSLELRARTVVEMLLHRHYLLQGMIDSKVMLQAVRSLCSE
 ESPEMTSLPPATLEAIFDADVKASCFSSFSNVWHLIALASVLQRNIYSIYPMRNLKIRPYFNVRIRPRR
 CDHVPSTLHIMWAGQPLTSHFFRHQYFAPVVGLEEEVEAEGAPGVAPALPALAPLSSPAKTLELLNREPL
 SYSHLCERYSVTKSTFYRWRQSQEHRQKVAARFSAKHFLQDSFHRGGVVPLQQFLQRFPEISRSTYYAW
 KHELLGSGTGPALPPREVLGMEELEKLPEEQVAEELECSALAVSSPGMVLQRAKLYLEHCISLNTLVP
 YRCFKRRFPGISRSTYYNWRKALRRNPSFKPAPALSAAGTPQLASVGEGAVIPWKSEAEAGNATGED
 PPAPGELLPLRMPLSRWQRRLRRAARRQVLSGHLPCFRFLRYPSLSPSAFVWVWVSLARGWPRGLSKLQV
 PVPTLGKGGQEAEEKQEKEAGRDVTAVMAPPVGASSEDVEGGSREGALQEGATAQGQPHSGPLLSQPVV
 AAAGGRDGRMLVMDMIATTKFKAQAKLFLQKRFQSKSFPSYKEFSALFPLTARSTYYMWKRALYDGLTLV
 DG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6689_d08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_018228


ORF Size: 2106 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_018228.3</u>
RefSeq Size:	3464 bp
RefSeq ORF:	2109 bp
Locus ID:	55237
UniProt ID:	<u>Q9H8Y1</u>
Cytogenetics:	14q24.3
MW:	78.3 kDa

Product images:



Circular map for RC208577



Western blot validation of overexpression lysate (Cat# [LY413208]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC208577 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).