

Product datasheet for **RC208605A1V**

Human CCR1 (NM_001295) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CCR1 (NM_001295) AAV Particle
Tag:	Myc-DDK
Symbol:	CCR1
Synonyms:	CD191; CKR-1; CKR1; CMKBR1; HM145; MIP1aR; SCYAR1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208605 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAACTCCAAACACCACAGAGGACTATGACACGACCACAGAGTTTGAATATGGGGATGCAACTCCGT
GCCAGAAGGTGAACGAGAGGGCCTTTGGGGCCCACTGCTGCCCTCTGTACTCCTTGGTATTTGTCAT
TGGCCTGGTTGGAAACATCCTGGTGGTCTGGTCTTGTGCAATACAAGAGGCTAAAAACATGACCAGC
ATCTACCTCCTGAACCTGGCCATTTCTGACCTGCTCTTCTGTTACGCTTCCCTTCTGGATCGACTACA
AGTTGAAGGATGACTGGGTTTTTGGTGATGCCATGTGTAAGATCCTCTGGGTTTTATTACACAGGCTT
GTACAGCGAGATCTTTTTCATCATCCTGCTGACGATTGACAGGTACCTGGCCATCGTCCACGCCGTGTTT
GCCTTGCGGGCACGGACCGTCATTTTGGTGTCATCACCAGCATCATCTTTGGGCCCTGGCCATCTTGG
CTTCCATGCCAGGCTTATACTTTTCCAAGACCAATGGGAATTCACCTACCACACCTGCAGCCTTCACTT
TCCTCACGAAAGCCTACGAGAGTGGAAGCTGTTTCAGGCTCTGAACTGAACCTCTTTGGGCTGGTATTG
CCTTTGTTGGTCATGATCATCTGCTACACAGGGATTATAAAGATTCTGCTAAGACGACCAATGAGAAGA
AATCCAAAGCTGTCCGTTTGATTTTGTATCATGATCATCTTTTCTCTTTTGGACCCCTACAATTT
GACTATACTTATTTCTGTTTCCAAGACTTCTGTTCACCCATGAGTGTGAGCAGAGCAGACATTTGGAC
CTGGCTGTGCAAGTGACGGAGGTGATCGCTACACGCACTGCTGTGTCAACCCAGTGATCTACGCCTTCG
TTGGTGAGAGGTTCCGGAAGTACCTGCGCAGTTGTTCCACAGGCGTGTGGCTGTGCACCTGGTTAAATG
GCTCCCCTTCTCTCCGTGGACAGGCTGGAGAGGGTCAGCTCCACATCTCCCTCCACAGGGGAGCATGAA
CTCTCTGCTGGGTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208605 protein sequence Red=Cloning site Green=Tags(s) METPNTTEDYDTTTEFDYGDATPCQKVNERAFGAQLLPPLYSLVFVIGLVGNILVVLVLVQYKRLKNMITS IYLLNLAISDLLFLFTLPFWIDYKLKDDWVFGDAMCKILSGFYTGLYSEIFFIILLTIDRYLAIVHAVF ALRARTVTFGVITSIIWALAILASMPGLYFSKTQWEFTHHTCSLHFPHESLREWKLQALKLNLFGLVL PLLVMIICTGTGIKILLRRPNEKKS KAVRLIFVIMIIFFLFWTPYNLTILISVFQDFLF THECEQSRHLD LAVQVTEVIAYTHCCVNPVIYAFVGERFRKYLRQLFHRRVAVHLVKWLPFLSVDRLERVSTSPSTGEHE LSAGF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001295
ORF Size:	1065 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001295.2
RefSeq Size:	2690 bp
RefSeq ORF:	1068 bp
Locus ID:	1230
UniProt ID:	P32246
Cytogenetics:	3p21.31
MW:	41.2 kDa