

Product datasheet for **RC208625L4V**

SLC37A2 (NM_198277) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC37A2 (NM_198277) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC37A2
Synonyms:	pp11662
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198277
ORF Size:	1515 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208625).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198277.1 , NP_938018.1
RefSeq Size:	4196 bp
RefSeq ORF:	1518 bp
Locus ID:	219855
UniProt ID:	Q8TED4
Cytogenetics:	11q24.2
Protein Families:	Transmembrane
MW:	54.8 kDa



[View online »](#)

Gene Summary:

Inorganic phosphate and glucose-6-phosphate antiporter. May transport cytoplasmic glucose-6-phosphate into the lumen of the endoplasmic reticulum and translocate inorganic phosphate into the opposite direction. Independent of a luminal glucose-6-phosphatase. May not play a role in homeostatic regulation of blood glucose levels.[UniProtKB/Swiss-Prot Function]