

Product datasheet for **RC208633L4V**

ZNF488 (NM_153034) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF488 (NM_153034) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF488
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153034
ORF Size:	1020 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208633).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153034.2
RefSeq Size:	3513 bp
RefSeq ORF:	1023 bp
Locus ID:	118738
UniProt ID:	Q96MN9
Cytogenetics:	10q11.22
Domains:	zf-C2H2
Protein Families:	Transcription Factors
MW:	37 kDa



[View online »](#)

Gene Summary:

Transcriptional repressor. Plays a role in oligodendrocyte differentiation, together with OLIG2. Mediates Notch signaling-activated formation of oligodendrocyte precursors. Promotes differentiation of adult neural stem progenitor cells (NSPCs) into mature oligodendrocytes and contributes to remyelination following nerve injury.[UniProtKB/Swiss-Prot Function]