

Product datasheet for **RC208639L3V**

APCDD1 (NM_153000) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	APCDD1 (NM_153000) Human Tagged ORF Clone Lentiviral Particle
Symbol:	APCDD1
Synonyms:	B7323; DRAPC1; FP7019; HHS; HTS; HYPT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153000
ORF Size:	1542 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208639).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153000.3
RefSeq Size:	2579 bp
RefSeq ORF:	1545 bp
Locus ID:	147495
UniProt ID:	Q8J025
Cytogenetics:	18p11.22
Protein Families:	Transmembrane
MW:	58.8 kDa



[View online »](#)

Gene Summary:

This locus encodes an inhibitor of the Wnt signaling pathway. Mutations at this locus have been associated with hereditary hypotrichosis simplex. Increased expression of this gene may also be associated with colorectal carcinogenesis.[provided by RefSeq, Sep 2010]