

Product datasheet for **RC208659L4V**

MKLP1 (KIF23) (NM_004856) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MKLP1 (KIF23) (NM_004856) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KIF23
Synonyms:	CHO1; KNLS5; MKLP-1; MKLP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004856
ORF Size:	2568 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208659).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004856.5
RefSeq Size:	3415 bp
RefSeq ORF:	2571 bp
Locus ID:	9493
UniProt ID:	Q02241
Cytogenetics:	15q23
Domains:	kinesin
Protein Families:	Druggable Genome



[View online »](#)

MW: 98.1 kDa

Gene Summary: The protein encoded by this gene is a member of kinesin-like protein family. This family includes microtubule-dependent molecular motors that transport organelles within cells and move chromosomes during cell division. This protein has been shown to cross-bridge antiparallel microtubules and drive microtubule movement in vitro. Alternate splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2013]