

Product datasheet for **RC208674**

ApoER2 (LRP8) (NM_001018054) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ApoER2 (LRP8) (NM_001018054) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ApoER2
Synonyms:	APOER2; HSZ75190; LRP-8; MCI1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Protein Sequence: >RC208674 representing NM_001018054
Red=Cloning site Green=Tags(s)

MGLPEPGPLRLLALLLLLLLLLLLRLQHLAAAAADPLLGGQGPAKECEKDQFQCRNERCIPSVWRCDDED
DCLDHSDEDDCPKKTCA SDFTCDNGHCIHERWKCDGEEPCDGSDESEATCTKQVCPAEKLSGPTSHK
CVPASWRCDGEKDCEGGADEAGCATLCA PHEFQCGNRSCLAAVFCVCDGDDDCGDGSDERGCADPACGPRE
FRCGGDGGGACIPERWVCDRQFDCEDRSDEAAELCGRPGPGATSAPAACATASQFACRSGEVHLGWRC
GDRDCKDKSDEADCP LGTRGDEFQCGDGTCLVAIKHCNQEQDCPDGSDEAGCLQGLNECLHNNGGC
CTDLKIGFECTCPAGFQLLDQKTCGDI DECKDPDACSQICVNYKGYFKCECYPGYEMDLLTKNCKAAAGK
SPSLIFTNRHEVRRIDLVKRNY SRLIPMLKNVVALDVEVATNRIYWCDLSYRKIYSAYMDKASDPKEQEV
LIDEQLHSPEGLAVDWHKHIIYWTDSGNK TISVATVDGRRRTLFSRNLSEPRAI AVDPLRGFMYSDWG
DQAKIEKSGLNGVDRQTLVSDNIEWPNGITL DLLSQRLYWVDSKHLQLSSIDFSGGNRKTLSSTDFLSH
PFGI AVFEDKVFWDLENAIF SANRLNGLEISILAENLNPHDIVIFHELKQPRAPDACE LSVQPNGGC
EYLCLPAPQISSHSPKYTCACPD TMLGPMKRCYRAPQSTSTTTLASTMTRTPATTRAPGTTVHRSTY
QNHSTETPSLTA AVPSSVSVPRAPSI SPSTLSPATSNHSQHYANEDSKMGSTVTA AVIGIIVPIVVIALL
CMSGYLIWRNWKRNKTKSMNFDNPVYRKTTEEEDEDELHIGRTAQIGHVYPARVALSLEDDGLP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI


Cloning Scheme:




ACCN: NM_001018054

ORF Size:	2712 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001018054.3
RefSeq Size:	4351 bp
RefSeq ORF:	2715 bp
Locus ID:	7804
UniProt ID:	Q14114
Cytogenetics:	1p32.3
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
MW:	99.13 kDa
Gene Summary:	This gene encodes a member of the low density lipoprotein receptor (LDLR) family. Low density lipoprotein receptors are cell surface proteins that play roles in both signal transduction and receptor-mediated endocytosis of specific ligands for lysosomal degradation. The encoded protein plays a critical role in the migration of neurons during development by mediating Reelin signaling, and also functions as a receptor for the cholesterol transport protein apolipoprotein E. Expression of this gene may be a marker for major depressive disorder. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jun 2011]

Product images:



Circular map for RC208674



Western blot validation of overexpression lysate (Cat# [LY422708]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC208674 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).