

## Product datasheet for **RC208679**

### PCGF5 (NM\_032373) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PCGF5 (NM_032373) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PCGF5
Synonyms:	RNF159
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC208679 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTACCCAAAGGAAACACTTGGTGAAAGATTTTAATCCTTACATTACCTGCTATATCTGTAAGGGT  
ATCTGATCAAGCCAACAACAGTGACGGAATGCCTCCATACATTCTGTAAGACTTGTATTGTTGAGCACTT  
TGAAGATAGCAATGATTGCCAAGGTGTGGCAACCAAGTTCATGAGACAAATCCATTAGAAATGTTGAGG  
TTGGACAATACATTAGAGGAAATTATATTTAAGCTGGTCCCTGGACTACGAGAACAAGAACTTGAGCGTG  
AATCTGAATTTTGGAAAGAAAATAGCCTCAAGAAAATGGACAAGATGATACTTCAAAGCTGACAAACC  
GAAAGTAGATGAAGAAGGTGATGAAAATGAAGATGATAAAGATTATCACAGAAGTGACCCACAAATTGCT  
ATCTGTCTAGATTGTTTACGAAATAATGGGCAATCAGGGGACAATGTAGTAAAGGGTTTAAATGAAGAAAT  
TCATTGATGTTCTACACGTGTAAGTGTGGAACTATTAATAAATTTCTAAGTTTAAACTAAAATTCC  
AAGTTCCTATGAGTTGGATGTGCTGTGCAATGGTGAATATGGGGAAGGATCATACTATGGAATTCATC  
TACATGACAAGATGGCGACTAAGAGGCGAAAACCTTCGGTGTCTGAAGTCTCAGCTTCGCAAGTCTGCT  
CTCAGGATGGCCCTTTGTATCAGTCATACCCTATGGTACTGCAGTATCGACCAAGAATTGATTTTCGGT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC208679 protein sequence  
 Red=Cloning site Green=Tags(s)

MATQRKHLVKDFNPYITCYICKGYLIKPTTVTECLHTFCKTCIVQHFEDSNDPCRCGNQVHETNPLEMLR  
 LDNTLEEIIIFKLVPLREQELERESEFWKKNKPQENGQDDTSKADKPKVDEEGDENEDDKDYHRSDPQIA  
 ICLDCLRNNQSGDNVVKGLMKKFIRCSTRVTVTIKKFLSLKLLKLPSSYELDLVLCNGEIMGKDHTMEFI  
 YMTRWRLRGENFRCLNCSASQVCSQDGPLYQSYPMVLQYRPRIDFG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: [https://cdn.origene.com/chromatograms/mk6361\\_b10.zip](https://cdn.origene.com/chromatograms/mk6361_b10.zip)

Restriction Sites: SgfI-MluI

Cloning Scheme:




Plasmid Map:



ACCN: NM\_032373

<b>ORF Size:</b>	768 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_032373.5</a>
<b>RefSeq Size:</b>	7222 bp
<b>RefSeq ORF:</b>	771 bp
<b>Locus ID:</b>	84333
<b>UniProt ID:</b>	<a href="#">Q86SE9</a>
<b>Cytogenetics:</b>	10q23.32
<b>MW:</b>	29.7 kDa
<b>Gene Summary:</b>	Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development. PcG PRC1 complex acts via chromatin remodeling and modification of histones; it mediates monoubiquitination of histone H2A 'Lys-119', rendering chromatin heritably changed in its expressibility (PubMed:26151332). Within the PRC1-like complex, regulates RNF2 ubiquitin ligase activity (PubMed:26151332). Plays a redundant role with PCGF3 as part of a PRC1-like complex that mediates monoubiquitination of histone H2A 'Lys-119' on the X chromosome and is required for normal silencing of one copy of the X chromosome in XX females (By similarity).[UniProtKB/Swiss-Prot Function]

**Product images:**

Western blot validation of overexpression lysate (Cat# [LY410166]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC208679 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).