

## Product datasheet for **RC208679L3V**

### PCGF5 (NM\_032373) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | PCGF5 (NM_032373) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | PCGF5  |
| Synonyms:                 | RNF159   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_032373  |
| ORF Size:                 | 768 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC208679).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_032373.3</a>  |
| RefSeq Size:              | 7222 bp  |
| RefSeq ORF:               | 771 bp   |
| Locus ID:                 | 84333  |
| UniProt ID:               | <a href="#">Q86SE9</a>   |
| Cytogenetics:             | 10q23.32   |
| MW:                       | 29.7 kDa   |



[View online »](#)

**Gene Summary:**

Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development. PcG PRC1 complex acts via chromatin remodeling and modification of histones; it mediates monoubiquitination of histone H2A 'Lys-119', rendering chromatin heritably changed in its expressibility (PubMed:26151332). Within the PRC1-like complex, regulates RNF2 ubiquitin ligase activity (PubMed:26151332). Plays a redundant role with PCGF3 as part of a PRC1-like complex that mediates monoubiquitination of histone H2A 'Lys-119' on the X chromosome and is required for normal silencing of one copy of the X chromosome in XX females (By similarity).[UniProtKB/Swiss-Prot Function]