

Product datasheet for **RC208685A1V**

Human TMEM120A (NM_031925) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TMEM120A (NM_031925) AAV Particle
Tag:	Myc-DDK
Symbol:	TMEM120A
Synonyms:	NET29; TMPIT
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208685 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAGCCCCGCCCCGGGCCGCTGGGCGACTGCCTGCGGGACTGGGAGGATCTACAGCAGGACTTCC
AGAACATCCAGGAGACCCATCGGCTCTACCGCTGAAGCTGGAGGAGCTGACCAAATTCAGAACAAATTG
CACCAGCTCCATCAGCGGCAGAGAAGCGGCTCCAGGAGCTGGCCCTCGCCCTGAAGAAATGCAAAACC
TCCCTCCCAGCAGAGGCCGAGGGGGCCGACAGGAGCTGGAGAACCAGATGAAAGAGCGCAAGGCCTCT
TCTTTGACATGGAGGCTATTTGCCTAAGAAGAATGGATTGTACCTGAGCCTGGTTCTGGGGAACGTCAA
CGTCACGCTCCTGAGCAAGCAGGCTAAGTTTGCCTACAAGGACGAGTATGAGAAGTTCAAGCTCTACCTC
ACCATCATCCTCATCCTCATCTCCTTCACTTGCCGCTTCTGCTCAACTCCAGGGTGACAGATGCTGCCT
TCAACTTCTGCTGGTCTGGTACTACTGCACCCTGACCATCCGGGAGAGCATCCTCATCAACAACGGCTC
CCGGATCAAAGGCTGGTGGGTGTTCCATCACTACGTGTCCGCTTCTGTGCGGAGTCATGCTGACGTGG
CCCGACGGTCTCATGTACCAGAAATCCGGAACCAATTCCTCTCCTTTTCCATGTACCAGAGCTTCGTGC
AGTTTCTCCAGTACTACTACCAGAGCGGCTGCCTCTACCGCTGCGGGCGCTGGGCGAGCGGCACACCAT
GGACCTCACTGTGGAGGGCTTCCAGTCTGGATGTGGCGGGGCTCACCTTCTGCTGCCTTTTCTTTTC
TTTGACACTTCTGGCAGCTTTTAACGCGCTGACGTTGTTCAACCTGGCCAGGACCCTCAGTGAAGG
AGTGGCAGGTGCTTATGTGCGGCTTTCCCTTCTCCTCTTTCTCGGCAATTTCTTACCACCCTGAG
GGTTGTGACACCACAAGTTTCACAGTCAGCGGCACGGGAGCAAGAAGGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208685 protein sequence Red=Cloning site Green=Tags(s) MQPPPPGPLGDCLRDWEDLQQDFQNIQETHRLYRLKLEELTKLQNNCTSSITRQKKRLQELALALKKCKP SLPAEAEAGAAQELNQMKERQGLFFDMEAYLPKKNGLYLSLVLGNNVTLLSKQAKFAYKDEYEKFKLYL TIILILISFTCRFLLSRVTDAAFNFLVWYYCTLTIREASILNNGSRIKGWVVFHHYVSAFLSGVMLTW PDGLMYQKFRNQFLSFSMYQSFVQFLQYYYQSGCLYRLRALGERHTMDLTVEGFQSWMWRLTFLLPFLF FGHFWQLFNALTFLNLAQDPQCKEQVLMCGFPFLLLFLGNFFTTLRVVHHKFHSQRHGSKKD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_031925
ORF Size:	1029 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_031925.1</u>
RefSeq Size:	1497 bp
RefSeq ORF:	1032 bp
Locus ID:	83862
UniProt ID:	<u>Q9BXJ8</u>
Cytogenetics:	7q11.23
MW:	40.6 kDa