

Product datasheet for RC208720L1

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GABPB2 (GABPB1) (NM_005254) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: GABPB2 (GABPB1) (NM_005254) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: GABPB2

Synonyms: BABPB2; E4TF1; E4TF1-47; E4TF1-53; E4TF1B; GABPB; GABPB-1; GABPB2; NRF2B1; NRF2B2

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC208720).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_005254

ORF Size: 1185 bp





GABPB2 (GABPB1) (NM_005254) Human Tagged Lenti ORF Clone - RC208720L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 005254.4</u>

 RefSeq Size:
 3069 bp

 RefSeq ORF:
 1188 bp

 Locus ID:
 2553

 UniProt ID:
 006547

Cytogenetics: 15q21.2

Domains: ANK

Protein Families: Transcription Factors

MW: 42.5 kDa


Gene Summary: This gene encodes the GA-binding protein transcription factor, beta subunit. This protein

forms a tetrameric complex with the alpha subunit, and stimulates transcription of target genes. The encoded protein may be involved in activation of cytochrome oxidase expression and nuclear control of mitochondrial function. The crystal structure of a similar protein in mouse has been resolved as a ternary protein complex. Multiple transcript variants encoding


distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC208720L1



Double digestion of RC208720L1 using Sgfl and Mlul $\,$