

Product datasheet for **RC208735L1V**

Calcium independent Phospholipase A2 (PLA2G6) (NM_001004426) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Calcium independent Phospholipase A2 (PLA2G6) (NM_001004426) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Calcium independent Phospholipase A2
Synonyms:	CaI-PLA2; GVI; INAD1; iPLA2; IPLA2-VIA; iPLA2beta; NBIA2; NBIA2A; NBIA2B; PARK14; PLA2; PNPLA9
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001004426
ORF Size:	2256 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208735).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001004426.1
RefSeq Size:	3077 bp
RefSeq ORF:	2259 bp
Locus ID:	8398
UniProt ID:	O60733
Cytogenetics:	22q13.1



[View online »](#)

Protein Pathways: alpha-Linolenic acid metabolism, Arachidonic acid metabolism, Ether lipid metabolism, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Glycerophospholipid metabolism, GnRH signaling pathway, Linoleic acid metabolism, Long-term depression, MAPK signaling pathway, Metabolic pathways, Vascular smooth muscle contraction, VEGF signaling pathway

MW: 84.1 kDa

Gene Summary: The protein encoded by this gene is an A2 phospholipase, a class of enzyme that catalyzes the release of fatty acids from phospholipids. The encoded protein may play a role in phospholipid remodelling, arachidonic acid release, leukotriene and prostaglandin synthesis, fas-mediated apoptosis, and transmembrane ion flux in glucose-stimulated B-cells. Several transcript variants encoding multiple isoforms have been described, but the full-length nature of only three of them have been determined to date. [provided by RefSeq, Dec 2010]

Product images:



[RC208735L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC208735L1V particle to overexpress human PLA2G6-Myc-DDK fusion protein.