

Product datasheet for **RC208780A1V**

Human UBE2L3 (NM_003347) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human UBE2L3 (NM_003347) AAV Particle
Tag: Myc-DDK
Symbol: UBE2L3
Synonyms: E2-F1; L-UBC; UBCH7; UbcM4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208780 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGCGGCCAGCAGGAGGCTGATGAAGGAGCTTGAAGAAATCCGCAAATGTGGGATGAAAACTTCCGTA
ACATCCAGGTTGATGAAGCTAATTTATTGACTTGGCAAGGGCTTATTGTTCTGACAACCTCCATATGA
TAAGGGAGCCTTCAGAATCGAAATCACTTTCCAGCAGAGTACCCATTCAAACCACCGAAGATCACATTT
AAAACAAAGATCTATCACCCAAACATCGACGAAAAGGGCAGGTCTGTCTGCCAGTAATTAGTGCCGAAA
ACTGGAAGCCAGCAACCAAAACCGACCAAGTAATCCAGTCCCTCATAGCACTGGTGAATGACCCCGAGCC
TGAGCACCCGCTTCGGGCTGACCTAGCTGAAGAATACTCTAAGGACCGTAAAAAATTCTGTAGAATGCT
GAAGAGTTTACAAAGAAATATGGGGAAAAGCGACCTGTGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208780 protein sequence
Red=Cloning site Green=Tags(s)

MAASRRLMKELEEIRKCGMKNFRNIQVDEANLLTWQGLIVPDNPPYDKGAFRIEINFPAEYPFKPPKITF
KTKIYHPNIDEKGQVCLPVISAENWKPATKTDQVIQSLIALVNDPQPEHPLRADLAEYSKDRKKFCKNA
EEFTKKYGEKRPVD

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_003347
ORF Size:	462 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003347.2</u>
RefSeq Size:	3028 bp
RefSeq ORF:	465 bp
Locus ID:	7332
UniProt ID:	<u>P68036</u>
Cytogenetics:	22q11.21
MW:	17.9 kDa