

## Product datasheet for **RC208795L4V**

### MTA3 (NM\_020744) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MTA3 (NM_020744) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MTA3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020744
ORF Size:	1545 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208795).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020744.2</a>
RefSeq Size:	2218 bp
RefSeq ORF:	1548 bp
Locus ID:	57504
UniProt ID:	<a href="#">Q9BTC8</a>
Cytogenetics:	2p21
Protein Families:	Transcription Factors
MW:	58.8 kDa



[View online »](#)

**Gene Summary:**

Plays a role in maintenance of the normal epithelial architecture through the repression of SNAI1 transcription in a histone deacetylase-dependent manner, and thus the regulation of E-cadherin levels. Contributes to transcriptional repression by BCL6.[UniProtKB/Swiss-Prot Function]