

## Product datasheet for **RC208815L4V**

### **NDNL2 (NSMCE3) (NM\_138704) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	NDNL2 (NSMCE3) (NM_138704) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NDNL2
Synonyms:	HCA4; LIC5; MAGEG1; MAGEL3; NDNL2; NSE3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138704
ORF Size:	912 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208815).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138704.3</a>
RefSeq Size:	1748 bp
RefSeq ORF:	915 bp
Locus ID:	56160
UniProt ID:	<a href="#">Q96MG7</a>
Cytogenetics:	15q13.1
Domains:	MAGE
MW:	34.3 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene is part of the SMC5-6 chromatin reorganizing complex and is a member of the MAGE superfamily. This is an intronless gene. [provided by RefSeq, May 2011]