

Product datasheet for **RC208821A1V**

Human DCUN1D4 (NM_015115) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human DCUN1D4 (NM_015115) AAV Particle
Tag: Myc-DDK
Symbol: DCUN1D4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208821 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCACTCGGATGCCGCCGTGTCAATTTTCAGCTGAACCTCTCATCTCTCAACACTGGCAAATATTCATA
AGATCTACCACACCCCTTAATAAGCTGAACCTAACAGAAGACATTGGCCAAGACGATCACCAAACAGGAAG
TCTGCGGTCTTGCACTTCTTCAGACTGCTTTAATAAAGTGATGCCACCAAGGAAAAAGAGAAGACCTGCC
TCTGGAGATGATTTATCTGCCAAGAAAAGTAGACATGATAGCATGTATAGAAAAATGATTTCGACTAGAA
TAAAGACTGAAGAAGAAGCCTTTTCAAGTAAAAGGTGCTTGAATGGTTCTATGAATATGCAGGAAGTGA
TGATGTTGTAGGCCCTGAAGGCATGGAGAAATTTTGTGAAGACATTGGTGTGAACCAAGAAAACGTAGTT
ATGCTTGTCTAGCTTGAAATTTGGATGCACAAAACATGGGTTATTTTACTCTACAGGAGTGGTAAAAAG
GAATGACTTCTCTCCAATGTGATACAACAGAAAACTCAGAAATACTTTGGATTACTTAAGATCATTCTT
AAATGATTCTACAACTTTAACTTATTTACAGATATGCGTTTGACTTTGCACGGCAATCAAAATACAAA
GTTATTAATAAAGACAGTGGTGCAATGTCCTAGAGTTTAGCAGAACAAATTAATCTTGACCTCAGCAACT
ATGATGAAGATGGAGCATGGCCAGTTTTGTTGGACGAGTTTGTGGAGTGGTATAAAGACAAACAGATGTC
C

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208821 protein sequence Red=Cloning site Green=Tags(s) MHSDAAAVNFQLNSHLSTLANIHKIYHTLNKLNLTEDIGQDDHQTGSLRSCSSSDCFNKVMPPRKKRRPA SGDDLAKKSRHDSMYRKYDSTRIKTEEEAFSSKRCLEWFYFYAGTDDVVGPEGMEKFCEDIGVEPENNV MLVLAWKLDAQNMGYFTLQEWLKGMTSLQCDTTEKLRNTLDYLRSLNDSTNFKLIYRYAFDFARQSKYK VINKDQWCNVLEFSRTINLDLSNYDEDGAWPVLLDEFVEWYKDKQMS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015115
ORF Size:	771 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_015115.2
RefSeq Size:	4283 bp
RefSeq ORF:	774 bp
Locus ID:	23142
UniProt ID:	Q92564
Cytogenetics:	4q12
MW:	30 kDa