

Product datasheet for **RC208826A1V**

Human CSP (DNAJC5) (NM_025219) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CSP (DNAJC5) (NM_025219) AAV Particle
Tag: Myc-DDK
Symbol: CSP
Synonyms: CLN4; CLN4B; CSP; DNAJC5A; mir-941-2; mir-941-3; mir-941-4; mir-941-5; NCL
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208826 representing NM_025219
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGACCAGAGACAGCGCTCACTGTCTACCTCTGGGGAGTCATTGTACCACGTCCTTGGGTTGGACA
AGAACGCAACCTCAGATGACATTAAGTCTATCGGAAGCTTGCTTGAAATATCACCCGACAAGAA
CCCCGACAACCCGGAGGCCGCGGACAAGTTAAGGAGATCAACAACGCGCACGCCATCCTACGGACGCC
ACAAAAAGGAACATCTACGACAAGTACGGCTCGCTGGGTCTCTACGTGGCCGAGCAGTTTGGGAAGAGA
ACGTGAACACCTACTTCGTGCTGTCCAGCTGGTGGGCCAAGGCCCTGTTGTCTTCTGCGCCTCCTCAC
GTGCTGCTACTGCTGCTGTCTGTCTGTCTGCTGCTTCAACTGCTGCTGCGGGAAGTGTAAAGCCCAAGGCG
CCTGAAGGCGAGGAGACGGAGTTCTACGTGTCCCCGAGGATCTGGAGGCACAGCTGCAGTCTGACGAGA
GGGAGGCCACAGACACGCCGATCGTCATACAGCCGGCATCCGCCACCGAGACCACCCAGCTCACAGCCGA
CTCCACCCAGCTACCACACTGACGGTTCAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208826 representing NM_025219
Red=Cloning site Green=Tags(s)

MADQRQRSLSTSGESLYHVLGLDKNATSDDIKSYRKLALKYHPDKNPDPNEAADKFKEINNAHAILTDA
TKRNIYDKYGLGLYVAEQFGEENVNTYFVLSSWWAKALFVFCGLLTCCYCCCCLCCCFNCCGKCKPKA
PEGEETEFYVSPEDLEAQLQSDEREATDTPIVIQPASATETTQLTADSHPSYHTDGFN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_025219
ORF Size:	594 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_025219.1</u>
RefSeq Size:	3286 bp
RefSeq ORF:	597 bp
Locus ID:	80331
UniProt ID:	<u>Q9H3Z4</u>
Cytogenetics:	20q13.33
MW:	22 kDa